




PRODUCT CATALOGUE

PARAMOUNT

QUEST FOR BETTER HEALTHCARE



The Paramount Philosophy:

With a vision of being a global leader in the medical disposable industry, Paramount Surgimed believes in providing its customers a timely supply of quality products at a competitive price. The company came into existence in 1993 under the Grover Group of companies capitalizing on the rich and diverse experience of over 65 years. Now, the brand is symbolized with trust and quality in the medical industry. The company started with the manufacturing of Surgical Blades and is now a leading supplier of a huge range of disposable products.

Paramount exports to over 97 countries, has 40 distributors and 400 Dealers across India. Our State-of-the-Art production facility compares to any international standard and is accredited with ISO: 13485:2016, CE Certificate and US FDA Device Listing.

The Paramount Vision:

Paramount Surgimed and all its employees share the vision of their dynamic Managing Director, Mr. Shaily Grover to be a global leader in the Medical Industry. Mr. Grover has ensured that the Group's reputation amongst its associates is synonymous with top quality, updated technology, dedicated approach and customer friendliness.



Corporate Office – Okhla, New Delhi



Factory – Bhiwadi, Rajasthan



Factory – Bhiwadi, Rajasthan



Factory – Bawal, Haryana

Surgical Blades	02
Disposable Scalpels	04
Safety Scalpel	05
Universal Safety Scalpel	06
Gouge Blades & Handles	07
Chisel/Fine/Micro Surgical Blade	08
Stitch Cutters & Thumb Scalpel	09
Dermal Blade & Metal Handle	11
Biopsy Punch	12
Dermal Curette	14
Ophthalmic Knives	16
Myringotomy Knife	18
Skin Graft Blade	19
Blood Lancet, Twisted & Flyfix	20
Disposable Prep Razor	21
Face Mask 3 PLY/ Type IIR	22
Face Mask N95/ FFP2	23
Disposable Adult Diapers	24
Under Pad - AroKleen & PupPee	26
Bed Bath Wipes, Cleansing Wipes	27
Germ Protection Wipes, Baby Wipes	28
Baby Wipes with water, Pet Wipes	29
Lung Exerciser, Hot Water Bottle	30
Infrared Thermometer, Pulse Oximeter & Air Mattress	31
Certificates	32

Surgical Blades

Carbon Steel & Stainless Steel

Features

- Stainless Steel Surgical Blades and Carbon Steel Surgical Blades.
- Packed in peel open pouches with a VCI protective inner paper.
- Shelf life of 5 years - Sterile Product
- Stainless Steel (Grade F as per BSEN ISO 7153-1:2001)
- Carbon Steel (Grade- As per BS 2982:1992)
- Complies the Hardness requirement of BS 2982:1992 standard.
- Complies the fitment requirement of ISO 7740:1985 standard requirement.



FITMENT 3



FITMENT 4





Details:

Product Name: Surgical Blades
 Size of Inner Box (cm): 8.5 x 7.0 x 2.7
 Unit/Box: 100 Pcs
 Outer Boxes: 10 Boxes x 100 Units
 Export carton: 10 Boxes x 10 Outer
 Export Carton size (cm): 42.5 x 39 x 17
 Per 20 feet container: 825 Carton

Size	Carbon Steel	Stainless Steel	Fitment No
9	PCS9	PSS9	3
10	PCS10	PSS10	3
10A	PCS10A	PSS10A	3
10R	PCS10R	PSS10R	3
11	PCS11	PSS11	3
11K	PCS11K	PSS11K	3
11P	PCS11P	PSS11P	3
12	PCS12	PSS12	3
12D	PCS12D	PSS12D	3
13	PCS13	PSS13	3
14	PCS 14	PSS14	3
15	PCS15	PSS15	3
15B	PCS15B	PSS15B	3
15C	PCS15C	PSS15C	3
15D	PCS15D	PSS15D	3
15S	PCS15S	PSS15S	3
16	PCS16	PSS16	3
17	PCS17	PSS17	3
18	PCS18	PSS18	4
19	PCS19	PSS19	4
20	PCS20	PSS20	4
21	PCS21	PSS21	4
22	PCS22	PSS22	4
22A	PCS22A	PSS22A	4
23	PCS23	PSS23	4
24	PCS24	PSS24	4
24D	PCS24D	PSS24D	4
25	PCS25	PSS25	4
26	PCS26	PSS26	4
36	PCS36	PSS36	4
36D	PCS36D	PSS36D	4
40	PCS40	PSS40	-
40B	PCS40B	PSS40B	-
60	PCS60	PSS60	8
60B	PCS60B	PSS60B	8
70	PCS70	PSS70	8

Post Mortem Blades –Size 60

Features:

- Used in the fields of Pathology and Histopathology for autopsy
- It is used with handle No.8
- Available in both material carbon steel and stainless steel
- Box of 50
- Individual peel pouch



Post Mortem Blades –Size 40

Features:

- The PM40 range of blades are meant to be used with screw type heavy duty handles
- Heavy duty blade
- Box of 10




Disposable Scalpels

Stainless Steel blade

Features

- Pre-assembled stainless steel surgical blades mounted on ABS plastic handle.
- The device is used to prevent injury during fixing the blade on the metal handle.
- Fluted finger grips and the blades are safely protected by LDPE guard.
- Individually packaged in a peel pouch.
- Convenience of single-use.
- Size clearly marked on each scalpel blade on blister pouch.
- Each dispenser box contains 10 scalpels.
- Shelf life of 5 years - Sterile Product



Details:

Product Name: Disposable Scalpel
 Size of Inner Box (cm): 21 x 7 x 3
 Unit/Box: 10 Pcs
 Export carton: 100 Boxes
 Export Carton size (cm): 44 x 37 x 33
 Per 20 feet container: 370 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PDSS9	3	18	PDSS18	4
10	PDSS10	3	19	PDSS19	4
10A	PDSS10A	3	20	PDSS20	4
11	PDSS11	3	21	PDSS21	4
11K	PDSS11K	3	22	PDSS22	4
12	PDSS12	3	22A	PDSS22A	4
13	PDSS13	3	23	PDSS23	4
14	PDSS14	3	24	PDSS24	4
15	PDSS15	3	25	PDSS25	4
15B	PDSS15B	3	36	PDSS36	4
15C	PDSS15C	3			
16	PDSS16	3			
17	PDSS17	3			

PDSS



Safety Scalpel

Stainless Steel blade


Features

- It is a safe alternative to the traditional scalpel and fulfills the requirement to reduce sharp injuries.
- Ribbed on both side of the sheath which provide extra gripping for left handed as well as right handed surgeon.
- One-handed activation allowing the surgeon to focus on the surgical site.
- Audible clicks makes the scalpel more secure and protective.
- "Used & Un-Used" indications prevent contamination and transmission of blood borne disease.
- Permanent safety lock prevents accidental cuts.
- Latex Free.
- Measurement scale 0 to 5 cm is embossed on the handle.
- Ergonomically designed safety scalpel.
- ABS plastic handle and sheath gives extra strength to the product.
- Universal size handle that is ideal for kits, trays & packs.

Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	300 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
10	PNSF10	3	21	PNSF21	4
10A	PNSF10A	3	22	PNSF22	4
11	PNSF11	3	23	PNSF 23	4
11K	PNSF11K	3	24	PNSF24	4
13	PNSF13	3	25	PNSF25	4
15	PNSF15	3	36	PNSF36	4
15B	PNSF15B	3			
15C	PNSF15C	3	MSC	PNSFMSC	3
17	PNSF17	4			
18	PNSF18	4			
19	PNSF19	4			
20	PNSF20	4			



Universal Safety Scalpel

Stainless Steel blade

Features

- **Safety Mode:** When you retract back after fully extending the blade, it makes a click sound to denote it is in safety mode. Safety mode means that the scalpel is not fully extended. This feature avoids injuries when the blade is half way extended as it retracts back and does not pierce through. It has to be fully extended to use.
- **Safety Lock:** if you completely retract back it makes a loud click denoting it is in safety lock position. Safety lock position means that the blade cannot be reused once it is in lock position. This feature forces the doctor to dispose the product and minimizes the risk of transmissibility of blood-borne infectious agents.
- Compatible for a left handed as well as right handed surgeon.
- Safety feature in pre use and post use of product.
- One-handed activation allowing the clinician to focus on the surgical site.
- The retractable design protects the blade from exposure, minimizing the risk of accidental cuts.
- Universal size handle that fits all types of hand sizes, and ideal for kits, trays & packs.
- One handle size fits all blade sizes.

Details:

Product Name:	Universal Safety Scalpel
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	43 x 37 x 20.5
Per 20 feet container:	280 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PSF9	3	18	PSF18	4
10	PSF10	3	19	PSF19	4
10A	PSF10A	3	20	PSF20	4
11	PSF11	3	21	PSF21	4
11K	PSF11K	3	22	PSF22	4
13	PSF13	3	22A	PSF22A	4
14	PSF14	3	23	PSF 23	4
15	PSF15	3	24	PSF24	4
15B	PSF15B	3	25	PSF25	4
15C	PSF15C	3			
16	PSF16	3	MSC	PSFMSC	3
17	PSF17	3			



Gouge Blades

Carbon Steel & Stainless Steel

Description

These surgical blades are used for foot surgeries and podiatry purposes. Gouge Blades are used by various podiatrist for hand, feet and skin surgery. They are available in the following sizes: Nos. 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 and 18. Each blade is protected by VCI paper in a peel- open foil sachet. Gouge blades fit into handle no. 3 and special gouge handles.



Details:

Product Name: Gouge Blades No. 1 to 5
 Size of Inner Box (cm): 10 x 10 x 3.3
 Unit/Box: 100 Pcs
 Outer Boxes: 10 Boxes x 100 Units
 Export carton: 10 Boxes x 10 Outer
 Export Carton size (cm): 59.5 x 47 x 19.5
 Per 20 feet container: 425 Carton

Size	Product Code	Product Code	Size	Product Code	Product Code	Fitment No	Handle
	Carbon Steel	Stainless Steel					
1	PGCS1	PGSS1	1V	PGCS1V	PGSS1V	3	PGMH
2	PGCS2	PGSS2	2V	PGCS2V	PGSS2V	3	PGMH
3	PGCS3	PGSS3	3V	PGCS3V	PGSS3V	3	PGMH
4	PGCS4	PGSS4	1R	PGCS1R	PGSS1R	3	PGMH
5	PGCS5	PGSS5	2R	PGCS2R	PGSS2R	3	PGMH
6	PGCS6	PGSS6	3R	PGCS3R	PGSS3R	3	PGMH
8	PGCS8	PGSS8				3	PGMH
10	PGCS10	PGSS10				3	PGMH
12	PGCS12	PGSS12				3	PGMH
15	PGCS15	PGSS15				3	PGMH
18	PGCS18	PGSS18				3	PGMH

Details:

Product Name: Gouge Blades No. 6 to 18
 Size of Inner Box (cm): 13.8 x 10 x 3.3
 Unit/Box: 100 Pcs
 Outer Boxes: 10 Boxes x 100 Units
 Export carton: 10 Boxes x 10 Outer
 Export Carton size (cm): 57 x 39 x 30
 Per 20 feet container: 390 Carton



Chisel Micro Surgical Blades

Fine Stainless Steel

Description

FINE (CHISEL) Surgical Blades Stainless Steel are available in the following sizes: Nos. Each blade is protected by VCI paper in a peel-open foil sachet. Sterilized by Gamma Radiation of 2.5 M. Rads.

Details:

Product Name: Chisel Blades (Fine)
 Size of Inner Box (cm): 8.2 x 3.2 x 2.2
 Unit/Box: 25 Pcs
 Outer Boxes: 25 Boxes x 25 Units
 Export carton: 25 Boxes x 10 Outer
 Export Carton size (cm): 47.5 x 39 x 17
 Per 20 feet container: 825 Carton

Size	Product Code	Handle
61	PSM61	PSMH
61V	PSM61V	PSMH
62	PSM62	PSMH
62V	PSM62V	PSMH
63	PSM63	PSMH
64	PSM64	PSMH
65	PSM65	PSMH
66	PSM66	PSMH
67	PSM67	PSMH
68	PSM68	PSMH
69	PSM69	PSMH
69B	PSM69B	PSMH
90	PSM90	PSMH
91	PSM91	PSMH



Stitch Cutters

Stainless Steel & Carbon Steel

Description

Stitch Cutter is used for general suture removal. Appearing similar to a scalpel, this instrument requires no handle and is a simple and convenient, way to remove basic sutures. The slightly curved, flat, bull nose tip slides under the suture and when rotated, gently removes the suture from the site. Long stitch cutter is made of stainless steel and short stitch cutter is made of carbon steel. The mini stitch cutter is made of stainless steel can be used with handle No.3L, 7 & 9. Alternatively, it is also available on pre-mounted ABS disposable plastic handle.

Features

- 5 Year Shelf Life – Sterile Product
- Each blade is protected by VCI paper in a peel-open foil sachet.
- Convenience of single-use disposables.
- Sizes clearly marked on each blade.
- Each box contains 100 pcs

Size	Carbon Steel	Stainless Steel	Handle No.
LSC	–	PLSCSS	-
SSC	PSSCCS	–	-
MSC	–	PMSCSS	PCS3L

Details:

Product Name:	Mini Stitch Cutter
Size of Inner Box (cm):	8.5 x 7.0 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	42.5 x 39 x 17
Per 20 feet container:	825 Carton



Details:

Product Name:	Short Stitch Cutter
Size of Inner Box (cm.):	11 x 8 x 3.5
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm.):	46.5 x 40.0 x 25.5
Per 20 feet container:	550 Carton



Stitch Cutters

Stainless Steel & Carbon Steel

Details:

Product Name:	Long Stitch Cutter
Size of Inner Box (cm):	14.8 x 7 x 3.1
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	45 x 34 x 36
Per 20 feet container:	475 Carton



Stitch Cutter with Handle

Stainless Steel

Details:

Product Name:	Mini Stitch Cutter with Handle
Size of Inner Box (cm):	14.5 x 70 x 3.5
Unit/Box:	Bulk Pack of 50 Pcs
Outer Boxes:	Not Using
Export carton:	10 Bags of 50 Pcs x 10 Outer
Export Carton size (cm):	42 x 32 x 39
Per 20 feet container:	290 Carton



Thumb Scalpel

Stainless Steel


Description

Disposable thumb scalpel or pocket scalpel is used in place of a surgical blade. This reduces or even prevents the incidences of injuries that can occur during the process of fitting or removing contaminated blades from handles.

The disposable thumb scalpels have a protective cap/guard for the blade that can be used so the scalpel can be safely disposed in the operation theater/room disposal container.

Details:

Product Name:	Thumb Scalpel
Product Code:	PTSS
Size of Inner Box (cm):	15.5x 11.5 x 8
Unit/Box:	25 Pcs
Export carton:	20 Boxes x 25 Units
Export Carton size (cm):	60 x 33.5 x 21
Per 20 feet container:	700 Carton



Dermal Blade


Stainless Steel blade

Description

- The device is intended for cutaneous surgery giving smooth clean excisions with less tissue trauma. It is Used to remove surface protuberances, flat lesions, or deep lesions by saucer-shaped incision. Dermal blade is made of stainless-steel blade mounted on flexible plastic holder. The device is under category of surgically invasive device and for transient use. It is covering under surface device (nature of body contact) and contact with intact skin and mucosal membrane. The device is invasive device and for transient use.

Details:

Product Name:	Dermal Biopsy Blade
Product Code:	PDB
Size of Inner Box (cm):	15.5 x 11.5 x 11
Unit/Box:	20 Pcs
Export carton:	20 Boxes
Export Carton size (cm):	60 x 33.5 x 24
Per 20 feet container:	700 Carton



Metal Handle

Description

- Metal handles comes in different sizes and are used according to the applications. The most popular stainless steel handle sizes used are size #3 and size #4 handles.

Gouge Handle - PGMH



Chisel Handle - PSMH



Metal Handle Fit.3 - PSSH3



Metal Handle Fit.4 - PSSH4



Biopsy Punch

Stainless Steel blade

Description

Biopsy punch is a sharp rounded blade tip used to remove tissue to be tested during procedures. The instrument is sterile with a round plastic handle. The punch is available in a variety of tip sizes and is sold in both sterile single pack as well as in bulk quantities.


Finely ground seamless, stainless steel blade ensures safe, clean and easy sample extraction. The sleek ribbed handle clear view of the biopsy site, comfort and control.

Virgin grade ABS plastic provides strength and sturdiness.

Details:


Product Name:	Biopsy Punch
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	43 x 37 x 20.5
Per 20 feet container:	635 Carton

Size mm	Product Code	Product Code
1.0	BPS1.0	BPP1.0
1.5	BPS1.5	BPP1.5
2.0	BPS2.0	BPP2.0
2.5	BPS2.5	BPP2.5
3.0	BPS3.0	BPP3.0
3.5	BPS3.5	BPP3.5
4.0	BPS4.0	BPP4.0
5.0	BPS5.0	BPP5.0
6.0	BPS6.0	BPP6.0
7.0	BPS7.0	BPP7.0
8.0	BPS8.0	BPP8.0



BPS Features

- Ergonomically designed for precise skin excisions
- Unique, ribbed, easy-grip barrel and pressure assistance platform
- Sterile and ready to use.
- Wide range of sizes



BPP Features

- Premium quality punch
- Hard blister packaging
- Option to print custom brand on handle
- Pressure assistance platform and ribbed handle
- Micro grinding finish



Dermal Curette

Stainless Steel blade

Description

The curette has a firm place in the range of dermatological, diagnostic and therapeutic instruments. It is one of the most versatile instruments. For that reason, it is increasingly also used in other disciplines, such as plastic surgery.


The aim of a curette is the removal of a superficial skin changes with minimal damage to the surrounding tissue.

The curette is extremely sharp thus ensuring precise and gentle tissue removal. The outstanding cutting ability makes working with it much easier compared with conventional products.

Details:

Product Name:	Dermal Curette
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	475 Carton

Size mm	Product Code	Product Code	Product Code	Product Code
2.0	DCS2.0	DCP2.0	DCK2.0	DDC2.0
3.0	DCS3.0	DCP3.0	DCK3.0	DDC3.0
4.0	DCS4.0	DCP4.0	DCK4.0	DDC4.0
5.0	DCS5.0	DCP5.0	DCK5.0	DDC5.0
7.0	DCS7.0	DCP7.0	DCK7.0	DDC7.0




DCS Features

- Seamless stainless steel sharp cutting edge
- Easy dynamic gripping handle
- Sterile and ready to use
- Sharp side indicator
- Wide range of sizes



DCK Features

- Sturdy & Durable Stainless Steel
- Symmetrically Ground Beveled Edge
- Ultra Sharp Inside Grinding
- Sharp side indicator
- Ribbed on all sides
- Triangular Handle
- Sterile and ready to use
- Wide range of sizes




DCP Features

- Long 'Pen Type' Grip
- Branding Labeling area
- Ribbed handle for easy control
- Minimum damage to surrounding tissues
- Sharp side indicator
- Sterile and ready to use
- Wide range of sizes



DDC Features

- Seamless stainless steel sharp cutting edge
- One in metal with a sharp edge
- Sharp side indicator
- Covered by a protective cap
- Ultra Sharp Inside Grinding
- Sterile and ready to use
- Wide range of sizes





Ophthalmic Knives

Round-Wire Blades



Slit-Blades (Keratome) For Phaco Stab Incision

Product No. Description

- P-912301 2.3 MM Sharp Tip Straight
- P-912501 2.5 MM Sharp Tip Straight
- P-912601 2.65 MM Sharp Tip Straight
- P-912801 2.8 MM Sharp Tip Straight
- P-912901 3.0 MM Sharp Tip Straight
- P-913201 3.2 MM Sharp Tip Straight
- P-913501 3.5 MM Sharp Tip Straight
- P-915101 5.1 MM Sharp Tip Straight

Double Bevel

- P-912808 2.8 MM Sharp Tip Straight
- P-912908 3.0 MM Sharp Tip Straight
- P-913208 3.2 MM Sharp Tip Straight

Product No. Description

- P-912361 2.3 MM Sharp Tip Angled 45° 19 Gauge
- P-912561 2.5 MM Sharp Tip Angled 45° 19 Gauge
- P-912661 2.65 MM Sharp Tip Angled 45° 19 Gauge
- P-912861 2.8 MM Sharp Tip Angled 45° 19 Gauge
- P-912961 3.0 MM Sharp Tip Angled 45° 19 Gauge
- P-913261 3.2 MM Sharp Tip Angled 45° 19 Gauge
- P-913561 3.5 MM Sharp Tip Angled 45° 19 Gauge
- P-915061 5.1 MM Sharp Tip Angled 45° 19 Gauge

Double Bevel

- P-912868 2.8 MM Sharp Tip Angled 45° 19 Gauge
- P-912968 3.0 MM Sharp Tip Angled 45° 19 Gauge
- P-913268 3.2 MM Sharp Tip Angled 45° 19 Gauge



Slit-Blades (N-Larger) For I.O.L. Enlargins

Product No. Description

- P-914001 4.0 MM Blunt Tip Straight
- P-915201 5.2 MM Blunt Tip Straight
- P-915501 5.5 MM Blunt Tip Straight
- P-916001 6.0 MM Blunt Tip Straight
- P-916201 6.2 MM Blunt Tip Straight

Product No. Description

- P-914061 4.0 MM Blunt Tip Angled 45° 16 Gauge
- P-915561 5.5 MM Blunt Tip Angled 45° 15 Gauge
- P-916061 6.0 MM Blunt Tip Angled 45° 15 Gauge
- P-916261 6.2 MM Blunt Tip Angled 45° 15 Gauge



Crescent Blades For Groove Or Pocket Or Tunnel Incision

Product No. Description

- P-950001 Crescent Sharp Tip Straight

Double Bevel

- P-950004 Crescent Angled 45° 20 Gauge
- P-950005 Sharp Tip Straight

Product No. Description

- P-950002 Crescent Bevel Up Angled 45° 20 Gauge
- P-950003 Crescent Bevel Down Angled 45° 20 Gauge



Lance Tip Blades For Initial Incision

Product No. Description

- P-931501 15° Degree 24 Gauge
- P-933001 30° Degree 24 Gauge
- P-934501 45° Degree 24 Gauge

Round-Wire Blades



MVR Blades For Water Tight Self Sealing Incision

Product No. Description

- P-975559 MVR 24 Gauge
- P-975560 MVR 20 Gauge
- P-975561 MVR 19 Gauge

Product No. Description

- AC STILETIO BLADES
- For Blumenthal ar Side Part Incision
- P-985560 MVR 20 Gauge Angled 45°
- P-985561 MVR 19 Gauge Angled 45°



Miniature Blades Alternative To General Surgery For Precise Incision

Product No. Description

- SPOON BLADES
- For Groove or Pocket or Tunnel Incision

Product No. Description

- P-6821 3.0 MM Angled 60°
- P-6821E 3.0 MM Extended Neck Angled 60°



Miniature Blades Alternative To General Surgery For Precise Incision

Product No. Description

- SCLERAL BLADES
- For Limbal and Scleral Incision
- P-5700 Wide Tip Blade
- P-5710 Narrow Tip Blade



Details:

Product Name:	Ophthalmic Knives
Size of Inner Box (cm):	20 x 5 x 4
Unit/Box:	6 Pcs
Export carton:	25 Boxes x 2 Outer
Export Carton size (cm):	42.5 x 27.5 x 22.5
Per 20 feet container:	1000 Carton

Myringotomy Knife Blades

Description

Myringotomy is a surgical procedure of the eardrum or tympanic membrane. The procedure is performed by making a small incision with a myringotomy knife through the layers of tympanic membrane.

Features:

- Optimised micrometre-accurate precision grinding.
- Double edged allowing for cutting movements in both directions.
- Ultra-sharp cutting edges of superb durability.
- Consistently sharp, angled blade facilitates a clean cut.
- Spear tip, Narrow shaft.
- Lance tip, Downward cut, Narrow shaft .
- With and Without Handle

Details:

Product Name:	Myringotomy Knife Blade
Size of Inner Box (cm):	24.2 x 8 x 2.5
Unit/Box:	5 Pcs
Outer Boxes:	--
Export carton:	50 Boxes
Export Carton size (cm):	50.5 x 42 x 14
Per 20 feet container:	1000 Carton

Size	Product Code
Lance	PSSML
Spear	PSSMS



Skin Graft Blade

Carbon Steel

Description


Skin Graft Blade - Simplex, Duplex

Surgical Skin Graft Blades are 158 mm long, individually wrapped in VCI paper in single peel packets in packs of 10. These blades can be used with Watson, Braithwaite or Cobbett Knives.

Details:

Product Name:	Skin Graft Blades
Size of Inner Box (cm):	20.7 x 3.6 x 2
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	44 x 24.5 x 22.5
Per 20 feet container:	1100 Carton

Size	Carbon Steel	Stainless Steel
SG	PCSG	PSSG
SGS	PCSGS	PSSGS
SGD	PCSGD	PSSGD



Blood Lancet

Description

Supplied sterile and individually packed. Lancet comes in a pack of 200 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.



Details:

Product Name:	Blood Lancet
Product Code:	PSSBL
Size of Inner Box (cm.):	12.5 x 7 x 20
Unit/Box:	200
Export carton:	100 Boxes x 200 Units
Export Carton size (cm.):	44.5 x 38.5 x 14
Per 20 feet container:	1800 Cartons



Blood Lancet - Twisted

Description

Supplied sterile and individually packed. Lancet comes in a pack of 100 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.



Details:

Product Name:	Blood Lancet
Product Code:	PSSBLT
Size of Inner Box (cm.):	6.5 x 4.2 x 7.0
Unit/Box:	100
Export carton:	64 Boxes x 100 Units
Export Carton size (cm.):	30 x 27.5 x 17
Per 20 feet container:	1800 Cartons



Flyfix

Features

- Discard after single use.
- Convenient.
- Extra strip for strong adhesion
- Saves nursing time.
- Infusion plaster with wound pad.

Details:

Product Code:	PAFF-24
Size of Box (cm.):	12x11x10
Unit/Box:	50
Export carton:	30 Boxes x 50 Units
Export carton size (cm):	51 x 39 x 25



Disposable Prep Razor

Paraprep

Features

- Safe contoured handle with sure grip ridges provide precise fingertips control.
- Easy to handle lies in front of blades thus requiring only light gentle pull instead of pressure.
- Coated stainless steel blades shaves close without nick or irritation.
- Open design provides optimal view of shaving area.
- Compact and convenient packaging.
- Cost effective, no need for more than one razor per procedure.

Details:

Product Code:	PPDR20210320
Size of Poly Bag (Cm.):	9.5x5.7
Inner Box:	20x12x6.5
Unit/Box:	50
Export carton:	20 Box x 50 Units
Export Carton size:	41.5X27.5X35
Per 20 feet container:	650 Cartons



with Two Side Blades



For All Body Areas



Smooth Hair Removal




100% Shower Proof



Global Standard Manufacturing



Unique Skin Protection





Stainless Steel Blade

Details:

Product Code:	PPRD
Size of Poly Bag (Cm.):	8.5x15.5
Inner Box:	28.5x19X26
Unit/Box:	100
Export carton:	5 Box x 100 Units
Export Carton size:	59 X 49 X 28
Per 20 feet container:	650 Cartons



Face Mask - 3 Ply/ Type IIR



3 PLY Mask - Non Surgical

3 PLY Mask - Surgical

- Smooth Fiber
Material
- SMS Fabric

- Meltblown
Fabric
- UV Cleaned

FEATURES

- Adjustable Nose
Clip
- Ultrasonic
Sealing
- One Size
Fits Most
- Seamless
Crimping
- Comfortable
& Breathable
- Quality
Material
- 3 Layer
Protection
- Anti Droplets

Details:

Product Name:	3 PLY/ Type IIR Mask - Non Surgical
Product Code	P03N
Size of Inner Box (cm.):	23 x 11 x 19.5
Unit/Box:	100
Export Carton size (cm.):	57 x 48 x 41
Per 20 feet container:	230 Carton

Details:

Product Name:	3 PLY/ Type IIR Mask - Surgical
Product Code	P03M
Size of Inner Box (cm.):	23 x 10.5 x 10
Unit/Box:	50
Export Carton size (cm.):	55 x 48 x 42
Per 20 feet container:	230 Carton

Face Mask - N95/FFP2



without Respirator
with earloop



UV Cleaned



without Respirator
with headloop



Filters Particulates



with Respirator
with earloop



BFE≥95%



PFE≥95%

FEATURES



Adjustable Noise Clip




Soft Elastic Ear Loops



High Quality Breathable Material and Valve



Seamless Flow of Exhaled Air



Details:

Product Name:	N95/FFP2 without Respirator
Product Code:	P95W
Size of Inner Box (cm.):	30.5 x 13.5 x 17.5
Unit/Box:	50
Export Carton size (cm.):	63.5 x 56 x 55
Per 20 feet container:	140 Carton

Details:

Product Name:	N95/FFP2 with Respirator
Product Code:	P95R
Size of Inner Box (cm.):	24.5 x 13.5 x 17.5
Unit/Box:	25
Export Carton size (cm.):	51 x 42 x 54.5
Per 20 feet container:	220 Carton

Disposable Adult Diapers For Urinary Incontinence Management

Features:

- Dry-Fast Core, Super absorbent.
- Moisture repellent.
- Security barriers.
- Waist Elastication and curved leg elastic.
- Superfit Band Wetness Indicator.
- Three absorbency levels and three sizes.
- Compact and convenient packaging.
- Also available in options of double cotton core & ADL.

Benefits:

- Faster absorption, drier surface.
- Better containment of faces.
- Reduces urine leakage.
- Improves comfort.
- Comfortable fit.
- Easier to fix and adjust tapes.
- A Guide for carers.
- Meets individual requirements with cost effective care.
- Easy to handle less storage space required.

Details:


Product Name:	Adult Diapers (Medium)
Product Code:	AKMP 04
Size of Poly Bag (cm.):	20 x 23.0 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	35.5x 46 x 64
Per 20 feet container:	240 Cartons

Details:

Product Name:	Adult Diapers (Large)
Product Code:	AKLP 06L
Size of Poly Bag (cm.):	22 x 24 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	48.5 x 35.5 x 68
Per 20 feet container:	220 Cartons

Details:

Product Name:	Adult Diapers (Extra Large)
Product Code:	AKXLP 07
Size of Poly Bag (cm.):	23.5 x 24 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	49 x 35.5 x 74
Per 20 feet container:	200 Cartons




Disposable Adult Diapers

Premium

-  Rapid Absorption Layer
-  Standing Leak Guards
-  Super Absorbent Anti-Bacterial Core
-  Odour Lock Technology
-  Additional Supporting Elastics

ADL (Acquisition Distribution Layer)



Details:

Product Code:	APM10
Size of Poly Bag (cm.):	15X18.5X18.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	51X47X36.5
Per 20 feet container:	290 Cartons

Economy

-  Rapid Absorption Layer
-  Standing Leak Guards
-  Super Absorbent Anti-Bacterial Core
-  Odour Lock Technology
-  Additional Supporting Elastics

Details:

Product Code:	APL11
Size of Poly Bag (cm.):	16X19X21
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	53X46.5X41
Per 20 feet container:	250 Cartons

Details:

Product Code:	APXL12
Size of Poly Bag (cm.):	20X16X21.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	54X48.5X42
Per 20 feet container:	230 Cartons



Under Pads

AroKleen

Description

- **Anti-Leakage Design:**
The product has extra width and length designed for better absorption.
- **Dry And Comfortable:**
It has large absorbing padding using the reticular embossed design to introduce the water diversion towards every side keep the surface comfortable.
- **Skin Friendly:**
Using the super-nonwoven as the top layer, it feels comfortable and prevents skin rashes.
- **International Quality:**
Strict controls of product material and production assuring safety and health.

Details:

Product Name:	Arokleen Baby Mat
Product Code:	ABM60x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Pkt:	10
Export carton:	15 Pkt x 10 Units
Export Carton size (cm.):	37 x 28 x 56
Per 20 feet container:	440 Cartons


Details:

Product Name:	Under Sheet
Product Code:	US90x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	41 x 27.5 x 55
Per 20 feet container:	410 Cartons

Details:

Product Name:	Puppy Pad
Product Code:	PP30P
Size of Poly Bag (cm.):	26.5 x 18 x 27
Unit/Pkt:	30
Export carton:	4 Pkt x 30 Units
Export Carton size (cm.):	60 x 36 x 29
Per 20 feet container:	290 Cartons

5 LAYERED LEAK PROOF TECHNOLOGY!



The Advantages of Using PupPee Pads

Our pads provide worry-free training for your puppy & peace of mind for you. They also can be used as an alternative to outdoor "relief" for your adult dog.

- Locks in wetness with quilted top layer
- Ensures maximum absorption power with 2 moisture-locking layers
- Protects floor & carpets with leak-proof backing
- Provides easy training with built-in attractant
- Protect your floors
- Easy cleanup & disposal
- Quilted for extra absorbency
- Less accidents means less disciplining
- Dog attractant reinforces training
- Durable, tear resistant



Bed Bath Wipes

AroKleen

Description

- Arokleen Bed Bath Wipes are Paraben Free, Alcohol Free, SLS Free, Mircowavable.
- The fast acting formula, with a two minute kill time loosens and suspends blood and bodily fluids, non-organic soil and other cellular debris for easy cleaning.

Details:

Product Code:	BP2202
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Pkt:	10
Export carton:	48 Pkt x 10 Units
Export Carton size(cm):	24x22x44
Per 20 feet container:	1100 Cartons

CLEANS, DISINFECTS & DEODORIZES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Skin Friendly pH



Aloe Vera



Enriched With Glycerin



USAGES



Cleansing



Refreshing



Moisturizing



Hospital Care

Cleansing & Refreshing Wipes

SPECIALLY FOR CLEANSING

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added Alcohol



Skin Friendly pH



Enriched With Vitamin C

Details:

Product Code:	BN2202
Size of Poly Bag (Cm.):	10X9.2X23.5
Unit/Pkt:	30
Export carton:	48 Pkt x 30 Units
Export Carton size(cm):	36 x 32 x 22
Per 20 feet container:	1000 Cartons



Germ Protection Wipes

AroKleen

Description

- Arokleen Germ Protection Wipes are more effective and cost less than sprays and towels, making them the right wipe to use in any location. Every area where staff, visitors, guests, customers and the general public convene, there are surfaces that need to be regularly cleaned and freed of dirt, grime, harmful germs and bacteria.

PROTECTION FROM GERMS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Skin Friendly
pH



Aloe Vera

USAGES



Hands



Refreshing



Restaurant



Car Interior



Device



Toys

Details:

Product Code:	PGP01
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Pkt:	48
Export carton:	36 Pkt x 48 Units
Export Carton size(cm):	36 x 33 x 23
Per 20 feet container:	920 Cartons



Baby Wipes

AroKleen

Description

- Arokleen Baby Wipes are large, durable, disposable clothes that are alcohol and latex free. Our Aloe Vera Extract, Vitamin E lotion-based formula that gently cleans and keeps babies' skin smooth and moisturized.

SPECIALLY FOR BABIES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added
Parabens



No Added
Alcohol



Extra Thick &
Extra Soft



Silicon Free



Enriched With
Vitamin E



pH
Balanced



Aloe Vera
Extract



Details:

Product Code:	BB23009
Size of Poly Bag (Cm.):	10X10.3X24.5
Unit/Pkt:	72
Export carton:	24 Pkt x 72Units
Export Carton size (cm):	37 x 32 x 25
Per 20 feet container:	850 Cartons



Baby Wipes with 99.9 % Water

SPECIALLY FOR KIDS
Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES

-  No Added Parabens
-  Aloe Vera
-  Silicon Free
-  pH Balanced
-  No Synthetic Preservatives

Details:

Product Code:	BWW01
Size of Poly Bag (Cm.):	10X10X24
Unit/Pkt:	72
Export carton:	24 Pkt x 72 Units
Export Carton size:	36 x 32 x 23.5
Per 20 feet container:	850 Cartons



Pet Wipes

SPECIALLY FOR CATS & DOGS
Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES

-  Powerful Odour
-  pH Balanced
-  Olive Oil
-  Enriched With Vitamin E
-  Cruelty Free
-  Vet & Groomer Recommended

Details:

Product Code:	PGP01
Size of Poly Bag (Cm.):	10X10X24
Unit/Pkt:	48
Export carton:	36 Pkt x 48 Units
Export Carton size:	38 x 32 x 22
Per 20 feet container:	950 Cartons



Lung Exerciser

Spirometer

Features:

- Three color balls for easier identification of inhale capacity
- Improve maximal exercise capacity
- Made from transparent material
- Compact, convenient and safe

Details:

Product Code:	PLE01
Size of Box (Cm.):	14.2x7x15.7
Unit/Box:	01
Export carton:	12 Bags x 01 Units



Hot Water Bottle



Joint Pain



Sports Pain



Arthritis Pain



Back Aches



Stomach Aches



Muscular Cramps



Rheumatic Pain



Muscle Pull



Regular Hot Water Bottle

Product Code: PHWB 05



BS Standard Hot Water Bottle

Product Code: HWBP 06

Infrared Forehead Thermometer

Features

- Non-contact Forehead
- LCD Backlight
- C/F Unit Switchable

Details:

Product Code:	PITM-10
Size of Box (Cm.):	10.5X4.65X16.5
Unit/Box:	01
Export carton:	12 Bags x 01 Units



16 Memories



1s Fast Reading



High Temperature Alarm



Pulse Oximeter

Features

- 2pcs standard AAA1.5V alkaline battery supporting 20 hours continuous operation.
- Applicable to Hospital, Clinic, Family, Oxygen Bar, Sporting center etc.
- True color high resolution OLED/TFT LCD/LED.
- Light in weight, simple in operation, easy to use & carry.
- Adjust the parameters in friendly menu.
- Soft anti-ambient light interference.
- Low power consumption.
- Optional Pouch.

Details:

Product Code:	PPOM-09
Size of Box (Cm.):	8.5x5x6
Unit/Box:	01
Export carton:	12 Bags x 01 Units



Air Mattress

Bedsore Prevention Cure Air Mattress

Features

- Anti decubitus alternating pressure air mattress Provide relief for bedsore
- Low noise, low vibration & low power Consumption.
- Recommended for bedridden patients.
- Medical grade PVC.
- Long working life.
- Large air output.

Details:

Product Code:	PACM-24
Size of Box(cm.):	26x16x20
Unit/Box:	01
Export carton:	06 Bags x 01 Units
Export carton Size:	50 x 28 x 42.5



16 Memories



High Temperature Alarm



1s Fast Reading



DRLM Home (mainMenu.htm)

Facility
 Products Listing

Annual Registration Successful

Facility: PARAMOUNT SURGIMED LTD., BHIWADI, Rajasthan, INDIA

You have successfully updated your registration and listing information for 2023.
 Your registration will be valid through Dec 31, 2023.
 Be sure to print this page for your records.
 The next registration renewal period is October 1 - December 31, 2023.
Registering your facility and listing devices does not, in any way, constitute FDA approval of your facility or devices.
 You may contact the FDA with any questions at regist@cdmr.fda.gov.
The Owner/Operator Number for this Registration is: 9012160.

Facility Information

Registration Number:
9680889
Initial Importer:
N
Facility Name:
PARAMOUNT SURGIMED LTD.
Legal Name:
Address:
A-106, RICO INDUSTRIAL AREA,
BHIWADI, Rajasthan, --, INDIA
DUNS Number:
Foreign Trade Zone:
N
Facility URL:
Other Business Trade Name(s):
Establishment located on a campus:

FDA

DNV

Certificate no.: 247943-2017-AQ-IND-NA-PS
 Place and date: Newark, 07 March 2022

Appendix to Certificate

Paramount Surgimed Ltd.
 Locations included in the certification are as follows:

Site Name	Site Address	Site Scope
Paramount Surgimed Ltd.	H.O.: 1, L.S.C., Okhla Industrial Area, Phase - II, New Delhi - 110 020, India	Design & Development, Sales & marketing of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fine Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel
Paramount Surgimed Ltd.	Site 2 : A-106, RICO Industrial Area, Bhowadi, District Alwar - 301 019, Rajasthan, India	Manufacture and Supply of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fine Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel
Paramount Surgimed Ltd.	Site 3 : Plot No. 33, Sector-14, Phase-2, IMT, Bawal - 123501, (Haryana), India	Manufacture and Supply of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fine Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel

Lack of Mitment of conditions as set out in the Certification Agreement may render this Certificate invalid.
 ACCREDITED UNIT: DNV Product Assurance AS, Verftveien 3, 1303 Høvik, Norway - TEL: +47 67 57 99 00, www.dnv.com

Page 2 of 2

ISO : 13485

Certificate of Registration

This is to certify that

Paramount Surgimed Ltd.
 H.O.: 1, L.S.C., Okhla Industrial Area, Phase II, New Delhi - 110 020, India.
 Works: A - 106, RICO Industrial Area, Bhiwadi, District Alwar - 301019 (Rajasthan), India.
 Works: Plot No. 33, Sector - 14, Phase - 2, IMT Bawal - 123501 (Haryana), India.

has been assessed by RIICL and found to comply with the requirements of

ISO 14001 : 2015
Environmental Management Systems

For the following activities:
 Design & Development, Manufacturing, Marketing and Supply of Medical Products like Adult Diapers, Disposable Undershorts, Baby Diapers, Sanitary Napkins, Under Pads, 2 Ply Face Masks, 3 Ply Face Masks, KN95 Face Masks, N95 Face Masks/ FFP2/FFP3, PPE Kit, Prep Razors, Baby Wipes, Bed Bath Wipes, Cleansing & Refreshing Wipes, Germ Protection Wipes, Disinfectant Wipes, Makeup Removal Wipes, Facial Cleansing Wipes, Dermal Blades and Section Needle.

This Certificate is Valid from 03/06/2023 Until: 02/06/2024

Date of Initial Certification: 03/06/2023
 1st Surveillance on or before: 02/05/2024
 1Ind Surveillance on or before: 02/05/2025
 Certification Valid Until: 02/06/2026

RIICL
 Director
 Royal Impact Certification Ltd.
 623, Tower-B, (Thum), Plot No. A-40, Sector-62, Noida 201301, India.
 www.riicl.com/india | info@riicl.com/india | Phone: +91 120 4113893

This certificate can be verified at www.riicl.com/india or www.riicl.com/research

UAF is member of International Accreditation Forum (IAF), the membership status can be verified at www.iaf.eu. This Certificate remains the property of Royal Impact Certification Limited. Must be returned on request or if certificate is withdrawn. Validity of this certificate is subject to successful surveillance audits as per dates given above

ISO : 4001

mna
 LABORATUVARLARI
 Nispetiye Bulvarı No: 2343

AB Tip İnceleme Sertifikası
EU Type-Examination Certificate

Belge No / Certificate No: 192062169
 Belgelendirme Tarihi - Bir Sınırlandırma Tarihi /
 Certification Date / Certificate Validity Date: 19.07.2023-03.07.2028
 Belge Geçerlilik Tarihi / Document Validity Period: 1.5 Yıl / 3 Years
 Firma Ünvanı ve Adresi /
 Company Name and Address: Paramount Surgimed Ltd
 A-106, RICO Industrial Area, Bhiwadi (Rajasthan)

Marka / Model / Brand / Model: 1. Aanklon
 Direktif / Direktive: 1. 2014/53 REGULATION
 Modül/Kategori / Module / Category: 1. B MODÜLÜ / KATEGORİ B)
 ABGÜRELE B / CATEGORY B)

Teknik Değerlendirme Raporu No /
 Technical Evaluation Report No: 1. MNA 92062169

Ürün Tipi / Product Type: 1. EN 149-2013 / A1 2009 Solunuma İlgili koruyucu ekipman - Filtresiz karp koruyucu amaçlı filtrelili yama maskesi / Respiratory protective devices - Filtering half masks to protect against particles

Ürünün Malzeme Bilgisi / Product Material Information: Aanklon modeli termal koruyucu, elastik kayış, burun klipsi ve filtre koruyucu kollektörün temel bileşenleri / Aanklon model products in one manufacturer's using fabric, elastic strap, nose clip filter layer

Yürürlükteki AB EN 149-2013 / 03.07.2023
 Karar Yürürlüğü / Approval: Şirket Müdürü / General manager

mna

MNA Laboratuvarları İnceleme Tarihi: 03.07.2023
 Adres: Nispetiye Bulvarı No: 2343 / Şişli / İstanbul / Türkiye / Tel: 0090 212 417 88 88 / Faks: 0090 212 417 19 19 / www.mna.com.tr

CE Mask



1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi-110 020, India

Tel +91-11-46436600, 46436601, 46436666, 26389812, 13

Fax +91-11-26389815, 41616555

sales@paramountblades.com

export2@paramountblades.com

www.paramountsurgimed.com

(Disclaimer- All measurements in weight & size are indicative for reference only.)

PSL/CTLG/01, Issue No. 01, Issue Date : 26/11/2020, Rev. No. 04, Rev. Date : 07/10/2024